


\* Determination soil reaction (pH)

Introduction:-

Soil reaction refers to the relationship between the concentration of hydrogen ions (H<sup>+</sup>) and hydroxyle ions (OH<sup>-</sup>) in the soil solution.



This ionization of water is controlled by ionization constant (K<sub>w</sub>)

$$K_w = \frac{(\text{H}^+) (\text{OH}^-)}{\text{HOH}} = 1 * 10^{-14}$$

The activity of HOH is equal 1 unit. In pure water the two ions are in equal concentration of 10<sup>-7</sup> mol/L solution

$$(\text{H}^+) (\text{OH}^-) = 1 * 10^{-14}$$


$$\text{And } (\text{H}^+) = (\text{OH}^-) = 1 * 10^{-7}$$

Sorenson 1909 defined pH as the **negative logarithm of the hydrogen activity**.

$$\text{pH} = - \log (\text{H}^+)$$

For example at pH= 5 this means that (H<sup>+</sup>) =10<sup>-5</sup> , and pH 6 means that H<sup>+</sup> =10<sup>-6</sup> , note that pH5 is 10 times more acidic than pH6

The following scale describe the pH scale ranges from 0 to 14



The pH of normal soils range between 3 and 9

### **Importance of soil reaction**

1- it affects the availability of nutrients in the soil. essential nutrients are most available to most plants at a pH between 6 to 7.5  
 at high pH, availability of phosphorus and most minor nutrients(except Boron and Molibidnium) decrease  
 at low pH, the solubility of the toxic heavy metals increase, and are absorbed by the plants

2- The CEC increase when pH increase, which also affect nutrients availability

### **Method of measurment of soil pH:**

the method of determination pH for soil are classified into :-

1- colorimetric method:-

- a) using litmos paper
- b) using indicators

the pH determination by this method involves the use of indicator dyes, which exhibit a change in color when the pH of a solution changes between certain limits.

The method is suitable for routine rapid pH determination, particularly under field conditions.

2- PH meter:-

### **Measuring the pH:-**

Soil suspension must be prepared by adding distilled water only

The ratio between Soil : Water is 1:1

Mix intermitently for 30 minutes then leave for 60 minutes

### **Apparatus**

pH meter with glass electrode

Glass beakers - 25 ml

### **Reagent**

Standard pH buffer solutions. Since soils of arid and semi-arid regions are alkaline in reaction, it is recommended to use, beside the pH 7 standard buffer solution, another standard solution in the pH range of 8.8 – 9.2 to calibrate the pH meter.

### **Procedure**

1. Standardize the pH meter using the standard pH buffers, with ample rinsing of the electrode with distilled water each time it is dipped into a buffer solution. Make sure to make adjustment for temperature correction according to instruction usually provided with the buffer.
2. Transfer a portion of the already prepared soil filtrate to fill about 3/4 of a 25 ml tall beaker.
3. Rinse the electrode carefully with distilled water and immerse into the soil filtrate.
4. Record the pH reading when the reading stabilizes.