

4- Complex Formation titration

Experiment no.11


Preparation and Standardization
of EDTA(di sodium salt)

Theory:

Chemical formula of Ethylene Diamine Tetra
Acetic acid di sodium salt is:



The Molar Mass is 372.25



Procedure:-

- 1- Prepare 0.01M EDTA(di sodium salt) 250mL
- 2- Prepare 0.01M Zinc Sulphate hepta hydrate (standard solution) 250mL
- 3-Take 10mL exactly Zinc Sulphate solution then add 40mL D.W.
- 4-Add 2mL pH10

5-Add small quantity of Erichrome BlackT

6-Titrate with EDTA solution in burette until the color of solution become blue.

7-Repeat 3,4,5and 6 twice

8-Calculate concentration of EDTA solution exactly.

*0.01M EDTA(di sodium salt) 250mL

Molar mass of EDTA(di sodium salt) = 372.25

Mass = 0.93 g

2- Prepare 0.01M Zinc Sulphate hepta hydrate

$\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ (standard solution) 250mL

Molar mass of $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O} = 287.44$

Mass = 0.71 g

At the endpoint

no. of millimoles of EDTA = no. of millimoles of
 $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$

$$M * V = M * V$$

$$M * Y = 0.01 * 10$$

$$M = ? \text{ Mol/L}$$