



## Experiment ( ): Qualitative Organic Analysis

( Identification of an Unknown: Aldehydes,  
Ketones & Carboxylic acids )

2022-2023

# Identification of an Unknown

Carboxylic acid

Ketone

Aldehyde




## CLASSIFICATION TESTS

which are simple chemical reactions that produce color changes or form precipitates, can be used to differentiate alcohols, aldehydes, and ketones and also to provide further structural information. Because color plays such an important role in this experiment.




The carbonyl group (C=O), found in aldehydes, ketones, is a very important functional group that is involved in several common reactions.


This particular functionality is unique because of the polarization (dipolar resonance) between the carbon-oxygen  $\pi$  bond.



carbonyl group



ketone



aldehyde

Due to the resonance of the C-O bond, there are a variety of reactions where the electrophilic carbon is attacked by nucleophiles (Lewis bases) and the oxygen reacts with corresponding electrophiles (Lewis acids).



Aldehydes commonly react faster than ketones (with the same reagent) due to less steric hindrance at the carbonyl group. Aldehydes are also more easily oxidized than ketones.

# 1- Identification of Aldehydes and Ketones

6

## 2,4-Dinitrophenylhydrazine (2,4-DNP) Test for Aldehydes and Ketones

Aldehydes and ketones undergo a condensation reaction with **2,4-dinitrophenylhydrazine** to produce yellow to orange precipitates as products. Alcohols do not undergo this reaction. Note the condensation reaction involves the removal of water with loss of the carbonyl oxygen.



## PROCEDURE:

-Clean and thoroughly rinse with distilled water four test tubes. Label a test tube for unknown (A, B, C, D).

1- Add 10 drops of each compound to its labeled test tube.

2- Add 5 drops reagent grade ethanol to each test tube.

3- Add 10 drops 2, 4- dinitrophenylhydrazine to each test tube.

Record your results. The presence of a yellow to orange-red precipitate is a positive test (aldehyde or ketone).

## 2- Tollen's Test


Tollens' test, also known as silver-mirror test, is a qualitative laboratory test used to distinguish between an aldehyde and a ketone. It exploits the fact that aldehydes are readily oxidized whereas ketones are not.

Tollens' test uses a reagent known as Tollens' reagent, which is a colorless, basic, aqueous solution containing silver ions coordinated to ammonia  $[\text{Ag}(\text{NH}_3)_2^+]$ . It is prepared using a two-step procedure.




**Step 1:** Aqueous silver nitrate is mixed with aqueous sodium hydroxide


9



**Step 2:** Aqueous ammonia is added drop-wise until the precipitated silver oxide completely dissolves.



Tollens' reagent oxidizes an aldehyde into the corresponding carboxylic acid.



Ketones are not oxidized by Tollens' reagent, so the treatment of a ketone with Tollens' reagent in a glass test tube does not result in a silver mirror.

### **PROCEDURE:**

First prepare the Tollen's reagent by:

- 1- Adding 5 mL of 5% silver nitrate to 2 freshly cleaned medium sized test tubes.
- 2-Add 3 drops of 10% sodium hydroxide to each test tube. Mix the solutions thoroughly. A brown-gray precipitate should form.
- 3-To each tube add 10% ammonium hydroxide solution drop-wise, shaking after the addition of each drop, until the precipitate just dissolves.


Then Add 2-3 drops of one of the unknowns to the 1<sup>st</sup> test tube and 2-3 drops of the other unknown to the 2<sup>nd</sup> one. Shake the tubes to mix them and allow them to stand without shaking for 10 minutes(heat if needed) . Look for the formation of the silver mirror.

### 3- Sodium Bicarbonate Test

11

Aldehydes and ketones are not acidic and will not react with sodium bicarbonate.

Carboxylic acid react with bicarbonate ion to produce carbon dioxide gas. This is another useful reaction, since it shows a visible change if the reaction occurs, you will see bubbles appear in the solution. This reaction is shown below.




## PROCEDURE

In this part, you will test your unknown add 1mL of 10% NaHCO<sub>3</sub> to test tube and then add 10 drops of unknown (Carboxylic acid , aldehyde or ketone).

While you are adding the substances to the sodium bicarbonate solution, watch for the formation of bubbles, which indicate a positive reaction (Carboxylic acid).

THANKS  
FOR YOUR ATTENTION

2,4-Dinitrophenylhydrazine (2,4-DNP or 2,4-DNPH) reacts readily with aldehydes and ketones via a condensation reaction (the lone pair of electrons on the terminal amino group in 2,4-DNPH makes it a strong nucleophile and the condensation starts by the nucleophilic 2,4-DNPH attacking the electrophilic carbonyl carbon) to produce the corresponding hydrazone. The hydrazone is usually a brightly colored yellow, orange or red compound, so the reaction is often used to test for the presence of an aldehyde or ketone



- 2,4-DNPH does not react with amides, esters or carboxylic acids. As shown below for the case of an ester, an extra resonance structure can be drawn for these 3 types of compounds as compared to a ketone. This extra resonance structure delocalizes some of the positive charge away from the carbonyl carbon onto the adjacent hetero-atom (oxygen in the case of the ester or carboxylic acid, nitrogen in the case of an amide). This makes the carbonyl carbon less electrophilic in these compounds, and consequently attack by the nucleophilic 2,4-DNPH is less favored.

